

Pancarpal arthrodesis in the dog

Professor John Innes BVSc PhD CertVR DSAS(orth) MRCVS
RCVS Specialist in small animal surgery (orthopaedics)



**VETERINARY
TISSUE BANK**

Indications

- 1 Palmar ligament rupture (traumatic and degenerative)
- 2 For radiocarpal bone fractures in breeds such as Boxers (Li and others 2000)
- 3 End-stage carpal arthritis (osteoarthritis, infective arthritis or erosive immune-mediated arthritis)
For non-reconstructable carpal fractures
- 4 As part of limb salvage for osteosarcoma of the distal radius and ulna

Patient preparation

The limb is clipped from just below the elbow joint to the level of the main carpal pad. An impervious barrier is taped in position over the toes. The patient is positioned in dorsal recumbency with the operative limb retracted caudally. Bilateral surgeries under the same anaesthetic are possible if required. The limb is free-draped and a sterile impervious drape wrapped over the foot. Sterilised cohesive dressing can then be tightly wrapped around the foot and progressed proximally to act as an Esmarch bandage; at the level of the proximal antebrachium, the dressing can be twisted in to a rope and used as a tourniquet.



Surgical approach

A skin incision is made through the cohesive dressing on the dorsal aspect of the carpus from the distal third of the antebrachium to the distal third of the third metacarpal bone. Care is taken to avoid the cephalic vein. The subcuticular tissues can then be stapled or sutured to the cohesive dressing on each side of the incision to isolate the wound.

Subcutaneous tissues are incised and the extensor carpi radialis tendon of insertion is identified on the medial aspect of the carpus. Lateral to this is the common tendon of the digital extensor muscle and this must be protected. Fascia between the ECR tendon and the CDE tendon is incised and the CDE tendon retracted laterally.

The joint capsule at all three joint levels is incised and the carpus flexed to open the joint space. Articular cartilage is removed using a combination of sharp dissection and curettage with manual or powered instruments. The subchondral plate of the radius is penetrated with 2-3 small drill holes to encourage mesenchymal cell migration.



Bone graft preparation

A suitable vial of Veterinary Tissue Bank freeze-dried canine cancellous chips (0.1-1mm) with a vial is chosen and rehydrated with sterile Hartmann's solution in a Gallipot. Graft volume is approximately 1cc for small dogs, 3cc for medium-sized dogs and 5cc for large/giant breeds.

Demineralized bone matrix (DBM) (Hoffer and others 2008) as available from Veterinary Tissue Bank may be mixed with the allograft chips to further enhance osteoinduction.